


Measure body temperature in seconds.  
Just by putting your face close to it.

## NINJA Mirror Thermometer

Non-contact emitter type Price: US\$145

Approach the mirror to a certain distance and align your forehead with the marking point. Only seconds without contact, by pressing a foot switch.

Entrance, walls or tables in hospitals, schools, train stations, stores, businesses or at home...



Prevents splash, no inhalation of splash.

## Vinyl 忍 Mask

A structure with filters and other components sandwiched between a rectangular piece of vinyl (PE), with multiple air vents around the nose and mouth.

About 90% of the surface of the mask is covered with PE vinyl, which prevents the diffusion and inhalation of splash, which cannot be prevented by ordinary paper masks.

Because using PE, it keeps the humidity inside the mask, moisturizes the throat and prevents the skin from drying out.

In addition, it uses less paper than a regular mask.

