

# NAOMi-CT Comparison of Scan Area & Object Size

201901

NEW

<b>Model</b>	<b>NAOMi-CT 001C</b>	<b>NAOMi-CT 002L</b>
<b>3D Panel</b>	Medium Panel	Large Panel
<b>Price</b>	US\$29,500	US\$39,500

\*Price does not include shipping fee.

SCAN AREA				
Scan Mode	Scan Area	Resolution	Scan Area	Resolution
<b>Offset Scan</b> (Wider Area)	φ151mm H82~63mm	0.205mm	φ240mm H175~125mm	0.24mm
<b>Normal Scan</b> (Normal)	φ 83mm H82~72mm	0.136mm	φ137mm H175~147mm	0.16mm
<b>Partial Scan</b> (High Resolution)	φ 56mm H41~37mm	0.068mm	φ68mm H49~45mm	0.083mm
Volume Comparison: <b>Approx. 5.16 times</b>  *With Offset Scan Area				

OBJECT SIZE (Size of the object that can be placed inside of the CT scanner)		
	φ250mm H230mm	φ300mm H320mm
Volume Comparison: <b>Approx. 2 times</b>		

SCANNER SIZE	

SCAN TIME	
About 1 min.	About 2 min.

\*Specifications and/or designs may be changed without a prior notice.