



Time bought back

My cost of developing, cost of chemical, processor maintenance and silver disposal went down to US\$300. On the contrary, my diagnostic ability becomes best. Now I can see clearer picture. Because of this I can treat patients well.

Built in cost of x-ray development is getting higher & higher. Also space for x-ray storage is huge. Carrying film & mailing cost is now zero.

Regular x-ray takes 3-5 min. to develop. NAOMI takes 30-60 sec. Immediate result. 3 min. saving X Number of x-ray = time bought back.


Patients were thrilled, more understanding about different management.



NOTE

No more development?

NAOMI uses sensor instead of film. Then the image is transferred to digital data and appear on the monitor shortly. Meaning? You can say goodbye to films, developer or even time spent on development. You don't even need CR, since NAOMI directly sends image to PC monitor.



PATIENTS WERE THRILLED