



Contact us

E-mail : [os@rfsystemlab.com](mailto:os@rfsystemlab.com)

# naomi-CT No X-ray Room Lead-shielded BOX

**Instant X-Ray Room**

No construction of special room

Ninja's CBCT fits anywhere!



# **NINJA Price** \$10,000

\*One price with no discount and no bargaining (straight game)  
\*Shipping is not included.


**No special lead room.  
Fits into any space in your clinic.**

Applicable with CT, Panoramic, and CT-Pano models.

BOX type CBCT contains lead shield, eliminating the needs for designated X-Ray room. You can install this anywhere, and enjoy your office space management.

## naomi-CT series **No X-ray Room Lead-shielded BOX Type**

width 1,220mm x Diphth 1,525mm x Height 1,900mm



Weight	250 kgs
Dimensions	W1,220 × D1,525 × H1,900 mm