



Contact us

E-mail : os@rfsystemlab.com

naomi-CT Dental Add-On

NINJA Price **\$1,900**

*One price with no discount and no bargaining (straight game)
*Shipping is not included.



Portable x-ray with flexible arm

The flexible arm enables you to position the x-ray to the patient, without asking him/her to move even after taking CT scan or panoramic x-rays.



Dental x-rays




Portable x-ray

naomi-CT series Dental Add-On Type

Width 1,220 mm
Depth 1,580 mm
Height 1,900 mm

*The arm bracket can be attached both on the left or right side.



X-ray	Tube voltage: 60 kV (Fixed)
	Tube current: 2 mA (Fixed)
Tube	Anode type: Fixed
	Target angle: 20 degrees
	Focal size: 0.8 mm
	Heat capacity 8.5 khu
Exposure time	0.01 - 1.6 sec (By 0.01 sec.)
Total iteration	1.6 mm Al or higher
Weight	2.0 kg
Diensions	W180 × D145 × H138 mm (Portable x-ray generator part)