



Contact us

E-mail : os@rfsystemlab.com

Dental Scanner




\$5,900

*One price with no discount and no bargaining (straight game)
*Shipping is not included.



■ Just like film, but easy and quick as this is "digital"

Only you need is to scan the IP (Imaging Plate) with a scanner, after taking x-rays just like films. The images will be displayed immediately. No troublesome with kids as IP is thin and small.



■ Specifications

IP size *1, 3, 4C are optional.	Size 0 (22 x 31 mm) 1 (24 x 40 mm)* 2 (31 x 41 mm) 3 (27 x 54 mm)* 4c (48 x 54 mm)*	Normal quality 343×484 pixels 375×625 pixels 484×640 pixels 421×843 pixels 750×843 pixels	High quality 628×885 pixels 685×1143 pixels 886×1171 pixels 771×1542 pixels 1370×1542 pixels
Pixel size	High quality: 35μm / Normal: 64μm		
Pixel pitch (35μm)	14.3 lp/mm		
Grayscale	Data acquisition: 8 bits		
Saving formats	DICOM/JPEG		
Weight	5.5kg		
Dimensions	120 × 318 × 265 mm		
Interface	High speed USB 2.0		
Operation system	Windows 7/8/8.1 and 10 *It does not work properly with some of Windows 10 series. *Not compatible with Apple Mac OS.		
Power	100-24V / 50-60Hz, 30W		
Attachments	IP Size 0 x2pcs, Size 2 x 4pcs / Disposable covers (100pcs each size) / IP case (1pc)		

*Product designs, specifications or features may be changed without a prior notice.