



Contact us E-mail : os@rfsystemlab.com

Wireless Intraoral Camera

Einstein Lumica



\$400

*One price with no discount and no bargaining (straight game)
*Shipping is not included.




Image available on a PC or TV, etc... with its designated receiver*.
(The receiver detects and display images from the wireless camera.)

■ Finger Tip Operation



■ Replaceable Battery



■ Specifications

Features	High-brightness white LED x 4 / Wireless / 64 pics on board / 1.36 megapixels
Transmission	Microwave transmitter
Processing method	Processing method
Charge	With the charging stand
Dimensions	φ25mm x 255mm
Weight	64.5g
Attachments	Battery pack, Charging stand, Disposable covers (100 sleeves), Lens protective sheet (5 sheets)

*Designated receiver (for PC or TV) is not included as an attachment. US\$200 each without shipping fee.

*Product designs, specifications or features may be changed without a prior notice.