



Contact us E-mail : [os@rfsystemlab.com](mailto:os@rfsystemlab.com)

# Portale X-ray ganerstor



# \$1,900

\*One price with no discount and no bargaining (straight game)  
\*Shipping is not included.



## Easy operations



## Specifications

X-ray	Tube voltage: 60 kV (Fixed)
	Tube current: 2 mA (Fixed)
Tube	Anode type: Fixed
	Target angle: 20 degrees
	Focal size: 0.8 mm
	Heat capacity 8.5 khu
Exposure time	0.01 - 1.6 sec (By 0.01 sec.)
Total iteration	1.6 mm Al or higher
Weight	2.0 kg
Diensions	W180 × D145 × H138 mm
Attachment	Charger

\*Product designs, specifications or features may be changed without a prior notice.