



Contact us

E-mail : os@rfsystemlab.com


# Slide CT-V

Large Size



# \$68,000

\*One price with no discount and no bargaining (straight game)  
\*Shipping is not included.




(Made in Japan)

■ No x-ray room required, it is a CT which can be installed with small space

**Anesthesia Free**  
No fears for dogs & cats  
with gently-covering air retention  
\*Attached with 3D Slide CT-V



CT Viewer  
忍者 NINJA




Simple process  
for scans



■ Dimensions

Slide CT-V

Width 1,030 × Depth 1,485 × Height 1,465 mm





Contact us

E-mail : [os@rfsystemlab.com](mailto:os@rfsystemlab.com)

## ■ Specifications

Specifications	Power	AC 100V 50/60Hz	
	Consumption power	2.0 kVA	
	X-ray generator	Tube voltage	60-90 kV
		Tube current	2-10 mA
		High voltage generator	DC rectification (Preignition method)
		Oil and air cooling	
		Focal spot	0.5 mm
		Total filtration	2.5 mm, Al or higher
	Emission switch	Dead man's switch	
	Dimensions	W1,030 mm × D1,485 mm × H1,465 mm	
	Weight	190 kgs	
Usage environment	Temperature	59~95°F / 15~35°C	
	Humidity	30~70%	
	Atmosphere pressure	920~1,060 hPa	
Partial mode	F.O.V	φ67 mm × H45~49 mm	
	Voxel size	0.083 mm	
Normal mode	F.O.V	φ141 mm × H147~176 mm	
	Voxel size	0.16 mm	
Wide mode	F.O.V	φ240 mm × H127~176 mm	
	Voxel size	0.24 mm	

\*Product designs, specifications or features may be changed without a prior notice.